

EXPERIMENTS WITH SINGLE TRAPPED YTTERBIUM IONS AT JPL

Nan Yu and Lute Maleki

Time and Frequency Sciences and Technology Group

Jet Propulsion Lab

California Institute of Technology


Pasadena, CA 91109

The work was carried out by the Jet Propulsion
Laboratory, California Institute of Technology, under a contract with the
National Aeronautics and Space Administration.



JPL

Ytterbium ion level scheme



369 nm cooling transition.



609nm 935nm, 638 nm
clearing transitions.

411nm, 467nm, 3.43um
clock transitions.



12 GHz microwave trans.

329 nm auxiliary shelving
transition.

Sketch of basic experimental setup



Paul-Straubel rf trap and single ion image



Single Yb⁺ ion fluorescence image





Photo of a Paul-Straubel trap

The trap is made of twisted Ta wire loop of 1.0 mm diameter.



The image resolution is limited by the photon collection optics (x4.5).

Resolved micro-motion sideband spectrum at 609 nm



609nm transition has a width of about 4MHz, conveniently resolving the 11MHz micromotion modulation sidebands. In the case, the micromotion amplitude is about 120nm.

D_{5/2} state lifetime measurement





The lifetime of the $D_{5/2}$ state can be determined using the quantum jump technique. In this technique, continuous fluorescence photons are detected when the ion is in the ground state. An excitation of the ion into the $D_{5/2}$ state will quench the fluorescence completely until it decays spontaneously back to the ground state and the fluorescence resumes.

The on or off state of the fluorescence signal indicates whether the ion is in the $S_{1/2}$ or $D_{5/2}$ state. The average fluorescence off-time(dark period) gives the lifetime of the $D_{5/2}$ state.


D_{3/2} state lifetime measurement

laser pulse sequence




The lifetime of D_{3/2} in Yb⁺ ion has been measured so far only with ion clouds in buffer gas, a condition not ideal for long lifetime measurement. It is not obvious how the measurement can be done with a single ion because the state is inside the transition cycle for the normal signal detection. *Our new technique takes advantage of the fact that the branching ratio of the P_{1/2} decaying to the ground state is large, about 200.*

measurement cycle:



D_{3/2} state lifetime measurement

decay histogram




On average, 207 photons at 369nm are scattered before the ion is repumped into D_{3/2} state again, yielding the P_{1/2} branching ratio to be 0.0483.

Motivations for individual ions in a high-finesse optical cavity

- Standing-wave laser cooling
- CQED in the strong coupling regime, structure and dynamics
- Real time probing of atom CM motion and measurement
- Quantum logic and information processing
- Quantum information distribution and networking
- Single ion laser/novel light source

Trapped individual ions in an optical cavity



Interaction Hamiltonian:

$$\hat{H}_s = \frac{\hbar\omega_A}{2}\hat{\sigma}^z + \hbar\omega_c\hat{a}^\dagger\hat{a} + i\hbar[g(\vec{r})\hat{a}^\dagger\hat{\sigma}^- - g^*(\vec{r})\hat{a}\hat{\sigma}^+].$$

$$g(\vec{r}) = \left(\frac{\mu^2\omega_c}{2\hbar\epsilon_0 V_m}\right)^{1/2}U(\vec{r}) = g_0U(\vec{r}).$$

The strong coupling condition:


$$g_0 > (\gamma, \kappa).$$

Experimental challenges:

- * protecting mirror coating,
- * avoid/circumvent dielectric charge-up,
- * reducing cavity volume,
- *

Initial exploratory system: special design points

- translatable trap/ion
- collimated atom beam
- collimated/rf-synchronized electron pulse
- separate cooling laser beam and cavity laser beam
- uv access for possible *insitu* surface discharge



Initial exploratory system: experimental goals


Feasibility demonstration:

- effects of atom beam collimation/contamination
- pulsed electron beam ion-loading/surface charge up problem
- ion translation capability, trap stability
- QED cavity locking/stabilization
- possible *insitu* surface discharge


Interesting physics to investigate:

- sw cooling
- cavity field mapping
- ion orbital size measurement
- laser transmission of occupied cavity
- QND atom state measurement through off resonance phase shift


Future systems: trap-cavity integration I




Future systems: trap-cavity integration II



linear trap ion loading and ion addressing



planar ring trap array



serial ring trap array