

Thursday Poster 083

LC/MS/MS as a Potential Method for Characterizing Bacterial Contamination During Spacecraft Assembly

Leticia Cornett¹; Cecilia Basic²; Wayne Schubert²,
Terry D. Lee¹

Beckman Research Institute, City of Hope, Duarte,
CA 91010¹

Jet Propulsion Laboratory, Pasadena, CA 91109²

Overview

Bacteria cells were processed by a lysis procedure followed by removal of cellular debris and digestion of the cytoplasmic components. LC/MS/MS analyses were performed using a custom built capillary HPLC system coupled to an ion trap mass spectrometer. MS/MS spectra were collected in a data dependent fashion and then matched to known protein sequences using the SEQUEST database search program.

Introduction

Exploration of outer bodies of interest to life's origins requires stringent measures to prevent contamination of these bodies with Earth life forms. Procedures are taken to vigorously clean spacecraft before launch. We are exploring MS to measure contamination during spacecraft assembly. MALDI and ESI-LC/MS create unique protein profiles from bacteria whole cells and cell lysates. Under ideal circumstances, these methods work well, but fail to yield consistent results on samples from different sources and handling techniques. Our more robust approach uses LC/MS/MS of the mixture of peptides obtained from the trypsin digestion of whole cell bacterial lysates

Methods


Freshly grown *B. subtilis* vegetative cell in 0.01% sodium azide were aliquoted into tubes and the solution was adjusted to 20% acetonitrile, 0.1% TFA. Cells were lysed by sonication (20 pulses 5 times at 1 minute intervals) on ice using a Branson sonifier 450. The lysate was centrifuged at 14,000 rpm for 30 minutes in the cold to remove cell debris then lyophilized. The proteins were denatured (urea), reduced (DTT) and alkylated (IAA) and dialyzed using a membrane with a 2K cutoff. An overnight trypsin digestion was performed. LC-MS-MS analysis was performed using a Finnigan LCQ, a custom built gradient loop HPLC and capillary (150u x 5cm) RP C18 column. All ion spectra were screened to eliminate poor MS/MS spectra using winnow program prior to database search.

Results

We are in the process of finding optimum conditions for the preparation of bacteria for MS analysis and the data presented are the combined results of 9 runs of varying conditions. Protein profiling did not produce reproducible results using bacteria prepared under different conditions and obtained from different sources. A trypsin digestion however allowed the positive identification of specific *Bacillus subtilis* proteins. Most of the proteins identified are expected to be present in high abundance in the cell. We found that the reduction/alkylation step was required to increase the number of hits by a factor of 2. Peak parking also increased the number of hits by a factor of 2.


Base Peak of sample. The analysis is fast and most peptides are concentrated in one region.

RT: 4.28 - 21.48




NL:
6.47E7
Base Peak


Example of a good MS/MS spectra that matched *B. subtilis* thioredoxin



MS/MS spectra that matched *B. subtilis* cold shock protein CSPB



MS/MS spectra that matched *B. subtilis* alanine dehydrogenase



protein	# times found	organism
owl f69756 f697	2	BACILLUS SUBTILIS
owl g69843 g698	4	BACILLUS SUBTILIS
owl g69985 g699	2	BACILLUS SUBTILIS
owl h69841 h698	2	BACILLUS SUBTILIS
owl O08327 GLGD_BACST	2	BACILLUS STEAROTHERMOPHILUS.
owl p00004 cyc_	7	contaminant
owl p02394 rl7_	8	BACILLUS SUBTILIS
owl p02662 cas1	8	contaminant
owl p02961 sasg	2	BACILLUS
owl p02968 fla_	5	BACILLUS SUBTILIS
owl p08877 pthp	9	BACILLUS SUBTILIS
owl p14949 thio	10	BACILLUS SUBTILIS
owl p26901 catv	6	BACILLUS SUBTILIS
owl p27876 tpis	4	BACILLUS SUBTILIS
owl p28599 ch10	3	BACILLUS SUBTILIS
owl p32081 cspb	5	BACILLUS SUBTILIS
owl p38494 rs1h	6	BACILLUS SUBTILIS
owl p41018 cspb	9	BACILLUS GLOBISPORUS
owl p49814 mdh_	8	BACILLUS SUBTILIS
owl p50727 fer_	3	BACILLUS SUBTILIS
owl p51777 cspd	4	BACILLUS SUBTILIS
owl p54542 yqje	7	BACILLUS SUBTILIS
owl p55179 pept	2	BACILLUS SUBTILIS
owl p80239 ahpc	10	BACILLUS SUBTILIS
owl p80698 tig_	4	BACILLUS SUBTILIS
owl p92505 cyc2	3	contaminant
owl p96707 ydgi	2	BACILLUS SUBTILIS
owl q08352 dha_	14	BACILLUS SUBTILIS
owl t01833 t018	2	Caenorhabditis elegans

definition

tellurium resistance protein homolog yceD - *Bacillus subtilis*

oligoendopeptidase homolog yjbG - *Bacillus subtilis*

3-hydroxybutyryl-CoA dehydratase homolog ysiB

conserved hypothetical protein yitW

GLYCOGEN BIOSYNTHESIS PROTEIN GLGD

contaminant

50S RIBOSOMAL PROTEIN L7/L12 (BL9) ('A' TYPE) (VEGETATIVE PROTEIN...

contaminant

SMALL, ACID-SOLUBLE SPORE PROTEIN GAMMA-TYPE (SASP)

FLAGELLIN.

PHOSPHOCARRIER PROTEIN HPR (HISTIDINE-CONTAINING PROTEIN).

THIOREDOXIN (TRX).

VEGETATIVE CATALASE

TRIOSEPHOSPHATE ISOMERASE (TIM).

10 KDA CHAPERONIN (PROTEIN CPN10) (PROTEIN GROES).

COLD SHOCK PROTEIN CSPB (MAJOR COLD SHOCK PROTEIN).

30S RIBOSOMAL PROTEIN S1 HOMOLOG.

COLD SHOCK PROTEIN CSPB

MALATE DEHYDROGENASE (VEGETATIVE PROTEIN 69) (VEG69)

FERREDOXIN

COLD SHOCK PROTEIN CSPD

HYPOTHETICAL 39.7 KD PROTEIN IN GLNQ-ANSR INTERGENIC REGION

PEPTIDASE T (AMINOTRIPEPTIDASE) (TRIPETIDASE).

ALKYL HYDROPEROXIDE REDUCTASE C22 PROTEIN (GENERAL STRESS PROTEIN 22).

TRIGGER FACTOR (TF) (VEGETATIVE PROTEIN 2) (VEG2).

contaminant

PUTATIVE NAD(P)H NITROREDUCTASE YDGI

ALANINE DEHYDROGENASE (STAGE V SPORULATION PROTEIN N)

wEST02554 Early embryo, *Stratagene* (cat. #937007) *Caenorhabditis elegans* cDNA clone CEESW23, mRNA

accession #

contaminant

contaminant

P02968

P08877

P14949

P26901

P27876

P28599

P32081

P38494

P41018

P49814

P50727

P51777

P54542

P55179

P80239

P80698

contaminant

P96707

Q08352

function

unknown; similar to tellurium resistance protein
unknown; similar to oligoendopeptidase
unknown; similar to 3-hydroxybutyryl-CoA dehydratase
unknown; similar to unknown proteins
required for glycogen biosynthesis
contaminant

contaminant

FORM THE FILAMENTS OF BACTERIAL FLAGELLA.
COMPONENT OF THE PHOSPHOENOLPYRUVATE-DEPENDENT SUGAR PHOSPHOTRANSFERASE SYST
REDOX REACTIONS & CATALYZES DITHIOL-DISULFIDE EXCHANGE REACTIONS
SERVES TO PROTECT CELLS FROM THE TOXIC EFFECTS OF HYDROGEN PEROXIDE.
PLAYS AN IMPORTANT ROLE IN SEVERAL METABOLIC PATHWAYS
BINDS TO CPN60 AND SUPPRESSES THE ATPASE ACTIVITY
CAN ACT AS TRANSCRIPTIONAL ACTIVATOR OF COLD SHOCK GENES
RIBOSOMAL PROTEINS
AFFECTS CELL VIABILITY AT LOW TEMPERATURES
[CATALYTIC ACTIVITY] $L\text{-MALATE} + NAD(+) = OXALOACETATE + NADH$
TRANSFER ELECTRONS IN A WIDE VARIETY OF METABOLIC REACTIONS.
[INDUCTION] IN RESPONSE TO LOW TEMPERATURE.
COULD BE A PEPTIDASE
HYDROLYZES A VARIETY OF TRIPEPTIDES CONTAINING N-TERMINAL METHIONINE, LEUCINE, OR PHE
DIRECTLY REDUCES ORGANIC HYDROPEROXIDES IN ITS REDUCED DITHIOL FORM.
INVOLVED IN PROTEIN EXPORT. ACTS AS A CHAPERONE BY MAINTAINING THE NEWLY SYNTHESIZED
contaminant
unknown; similar to NADH dehydrogenase
ASSIMILATION OF L-ALANINE AS AN ENERGY SOURCE THROUGH THE TCA CYCLE DURING

EM (PTS),

NYLALANINE

PROTEIN IN AN OPEN CONFORMATION

Example of a SEQUEST summary file grouped by consensus

SEQUENCE Summary - Microsoft Internet Explorer

File Edit View Favorites Tools Help

Back Forward Stop Refresh Home Search Favorites History Mail Print

Address http://151.152.4.251/cgi-bin/runsummary.pl?directory=icornetabridgestfilter&sort=consensus Go Links

#	rBP	Efile	e	dH	mK	AARS	SLS	Sp	PSP	Tons		Q	Sequence	Full to Top
<input checked="" type="checkbox"/> A owl P32081 CSPB_BACSU 30 6 3.2e6 (3,0,0,0,0,0) (4 35 37 41 43 44, x, x, x, x, x) COLD SHOCK PROTEIN CSPB (MAJOR COLD SHOCK PROTEIN). - BACILLUS...														
35	1.0e4	0579	2	0.4	2859.0	4.41	0.49	605	1	21/50	<input checked="" type="checkbox"/>	owl p41018 cspb +1	(K) GFGFIEVEGQDDVVFVHFSAIQGEGFK	
44	2.3e3	0636	3	-0.6	2859.9	3.65	0.38	1302	88	37/100	<input checked="" type="checkbox"/>	owl p41018 cspb +1	(K) GFGFIEVEGQDDVVFVHFSAIQGEGFK	
43	2.2e3	0630	2	-1.0	1878.9	2.88	0.41	1103	1	20/32	<input checked="" type="checkbox"/>	owl p32081 cspb	(K) TLEEGQAVSFEIVEGNR	
41	3.4e3	0618	2	1.4	2857.9	2.60	0.25	643	1	23/50	<input checked="" type="checkbox"/>	owl p41018 cspb +1	(K) GFGFIEVEGQDDVVFVHFSAIQGEGFK	
37	3.5e3	0600	2	-1.7	1879.6	2.42	0.24	369	265	15/32	<input checked="" type="checkbox"/>	owl p32081 cspb	(K) TLEEGQAVSFEIVEGNR	
4	3.1e6	0462	1	0.2	885.3	2.29	0.38	461	11	13/16	<input checked="" type="checkbox"/>	owl p32081 cspb +1	(R) GPQAANVTX	
<input checked="" type="checkbox"/> B owl Q08352 DHA_BACSU 30 3 2.6e6 (3,0,0,0,0,0) (15 19 39, x, x, x, x, x) ALANINE DEHYDROGENASE (EC 1.4.1.1) (STAGE V SPORULATION PROTEIN)...														
39	6.6e3	0612	3	-0.2	2961.9	3.40	0.31	1171	85	39/112	<input checked="" type="checkbox"/>	owl q08352 dha	(K) TLISNPVNIADAVAADLLICAVLIPGA	
19	2.1e6	0519	1	0.1	1138.6	2.19	0.23	171	2	15/22	<input checked="" type="checkbox"/>	owl q08352 dha	(K) GILLAGVPGVSR	
15	4.1e5	0504	1	0.3	1356.5	1.53	0.20	90	1	18/28	<input checked="" type="checkbox"/>	owl q08352 dha	(K) VTIIGGGVVGTNAAK	
<input checked="" type="checkbox"/> C owl P49814 MDH_BACSU 20 2 8.9e4 (2,0,0,0,0,0) (32 36, x, x, x, x, x) MALATE DEHYDROGENASE (EC 1.1.1.37) (VEGETATIVE PROTEIN 69) (VEG69)...														
32	7.7e4	0564	2	0.6	2558.8	4.92	0.43	1139	1	23/44	<input checked="" type="checkbox"/>	owl p49814 mdh	(K) YSPDSIIIVLTNTPVDAMTYAVYK	
36	1.2e4	0585	3	-0.2	3236.9	3.97	0.18	1528	7	42/124	<input checked="" type="checkbox"/>	owl p49814 mdh	(R) KGGGEIVNLLNGCSAYYAPAAASLTENVAILK	
<input checked="" type="checkbox"/> D owl P51777 CSPD_BACSU 20 2 1.5e5 (2,0,0,0,0,0) (8 40, x, x, x, x, x) COLD SHOCK PROTEIN CSPD. - BACILLUS SUBTILIS.														
40	3.0e3	0615	2	0.7	1921.2	3.10	0.16	955	23	17/32	<input checked="" type="checkbox"/>	owl p51777 cpspd	(K) SLEEGQEVSFEIVEGNR	
8	1.5e5	0474	1	0.2	899.3	1.73	0.30	399	369	11/16	<input checked="" type="checkbox"/>	owl p51777 cpspd	(R) GPQASNVTX	
<input checked="" type="checkbox"/> E owl P80239 AHPC_BACSU 20 2 7.8e4 (2,0,0,0,0,0) (1 27, x, x, x, x, x) ALKYL HYDROPEROXIDE REDUCTASE C22 PROTEIN (EC 1.6.4.-) (GENERAL...														
27	3.2e4	0543	2	0.9	2341.3	2.66	0.15	1010	1	24/44	<input checked="" type="checkbox"/>	owl p80239 ahpc	(R) GTFIIDPDGVIQTVEINAGGGR	
1	4.5e4	0453	1	0.2	945.3	2.24	0.29	507	38	10/14	<input checked="" type="checkbox"/>	owl p80239 ahpc	(K) GWHDSSEK	

Internet

Table of most frequently observed *B. subtilis* proteins and their function

Microsoft Excel - Table of common proteins & functions.xls		
File Edit View Insert Format Tools Data Window Help		
Arial 10 B I U E F 75%		
B26		
1	A	C
	protein name	protein function
2	tellurium resistance protein homolog yoeD	P63756 unknown
3	heavy metal-transporting ATPase homolog ykvV	P63669 unknown
4	oligoendopeptidase homolog yibG	G63643 unknown
5	probable enoyl-CoA hydratase (EC 4.2.1.17) ysfB	G63685 hydro-lyase
6	50S RIBOSOMAL PROTEIN L7/L12 (BL9) (A' TYPE) (vegetative protein 34)	P62384 protein synthesis
7	FLAGELLIN	P62368 FORMS THE FILAMENTS OF BACTERIAL FLAGELLA
8	phosphocarrier PROTEIN HPR (HISTIDINE-CONTAINING PROTEIN)	P63877 component of the phosphoenolpyruvate-dependent sugar phosphotransferase system
9	THIOREDOXIN (TRX)	P14343 VARIOUS REDOX reactions & CATALYZES DITHIOL-DISULFIDE EXCHANGE reactions
10	vegetative catalase	P26301 PROTECT CELLS FROM THE TOXIC EFFECTS OF hydrogen peroxide
11	TRIOSE-PHOSPHATE ISOMERASE	P27876 PLAYS AN IMPORTANT ROLE IN SEVERAL METABOLIC PATHWAYS
12	10 KDA CHAPERONIN (PROTEIN GROES)	P29559 SUPPRESSES THE ATPASE ACTIVITY OF CPN60
13	COLD SHOCK PROTEIN CSPB (MAJOR COLD SHOCK PROTEIN)	P32031 TRANSCRIPTIONAL ACTIVATOR OF COLD SHOCK GENES
14	36S RIBOSOMAL PROTEIN S16 HOMOLOG	P38434 BELONGS TO THE S16 FAMILY OF RIBOSOMAL PROTEINS
15	COLD SHOCK PROTEIN CSPB	P4108 AFFECTS CELL VIABILITY AT LOW TEMPERATURES
16	MALATE DEHYDROGENASE (VEGETATIVE P-PROTEIN 63)	P43814 [CATALYTIC ACTIVITY] L-MALATE + NAD(+) = OXALOACETATE + NADH
17	FERREDOXIN	P60727 TRANSFER ELECTRONS IN A WIDE VARIETY OF METABOLIC REACTIONS
18	COLD SHOCK PROTEIN CSPD	P51777 [INDUCTION] IN RESPONSE TO LOW TEMPERATURE
19	hypothetical 39.7 KD PROTEIN IN GLN-ANSP INTERGENIC REGION	P54542 COULD BE A PEPTIDASE
20	PEPTIDASE T (AMINOTRIPEPTIDASE) (TRIPEPTIDASE)	P55179 HYDROLYZES A VARIETY OF TRIPEPTIDES CONTAINING N-TERMINAL MET, LEU, OR PHE
21	ALKYL HYDROPEROXIDE REDUCTASE C22 PROTEIN (GENERAL STRESS PROTEIN 22)	P60239 DIRECTLY REDUCES ORGANIC HYDROPEROXIDES IN ITS REDUCED DITHIOL FORM
22	TRIGGER FACTOR (TF) (VEGETATIVE PROTEIN 2)	P60638 INVOLVED IN PROTEIN EXPORT. ACTS AS A CHAPERONE
23	PUTATIVE NAD(P)H NITROREDUCTASE YDG1	P36717 unknown
24	ALANINE DEHYDROGENASE (STAGE V SPOULATION PROTEIN N)	Q08352 KEY FACTOR IN THE ASSIMILATION OF L-ALANINE AS AN ENERGY SOURCE THROUGH THE TCA CYCLE DURING
25		
26		
27		
28		
29		
30		
31		
32		
33		
34		
35		
36		
37		
38		
39		
40		
41		
42		
43		
44		
45		
Sheet1 Sheet2 Sheet3		
Ready		

Conclusions

The described procedure depends on obtained fragment ion spectra of sufficient number of the most abundant and readily ionizable peptides present in the cytoplasm for that particular organism. The technique was found to work well for organisms whose genome has been sequenced. For *Bacillus subtilis*, whose genome is in the database, we obtained specific peptide matches to expected proteins such as chaperonins, metabolic enzymes and ribosomal proteins. While the complex bacterial digestion mixture did not yield the same set of results each time different aliquots or concentrations were analyzed, the majority of high scoring matches were assigned to *Bacillus subtilis* proteins.

We were also able to identify a few hypothetical proteins in the database previously identified by open reading frames. We are unsure at this time whether the presence of these proteins have been detected by any other means. As more bacterial genomes are sequenced, this procedure can be applied to a variety of applications including the analysis of bacteria in distinct microenvironments.

A number of high scoring hits for *C. elegans* was also observed and we later found that this organism is commonly used in the laboratory where the bacteria are being grown hence a possible contaminant for our preparations. The fact that we found this contaminant leads us to believe that our method is sensitive to the presence of other organisms and may be used to identify mixtures of heterogeneous organisms.

References

National Center for Biotechnology Information (NCBI) biomedical databases

Kunst et al (1997) "The complete genome sequence of the Gram-positive bacterium *Bacillus subtilis*"
Nature 390:249

Altschul et al (1997) "Gapped BLAST and PSI-BLAST: a new generation of protein database search programs", Nucleic Acids Res. 25:3389-3402.

Yates et al. (1998) " Method to compare collision-induced dissociation spectra of peptides: potential for library searching and subtractive analysis" Anal Chem 70:3557-65

Link et al (1999) "Direct analysis of protein complexes using mass spectrometry" Nat Biotechnol 1999 Jul;17(7):676-82

Krishnamurthy et al (1999) "Liquid chromatography/microspray mass spectrometry for bacterial investigations " *Rapid Commun Mass Spectrom* 1999;13(1):39-49

Moore et al. (2000) "Method for Screening Peptide Fragment Ion Mass Spectra Prior to Database Searching" JASMS 11:422-426

Acknowledgements

The authors wish to thank the following individuals:

Wayne Shubert for providing us the bacteria cultures.

Curtis Croker for constant modification of our software-controlled HPLC and the peak parking procedure which will be used for any further separations.

John Percy for his programming skills used for the creation of the MoWeD deconvolution software which was initially used in the protein profiling experiments.

Roger Moore for creation of the winnow program which was used as an effective screen for MS/MS spectra.

Acknowledgements

Curtis Croker for constant modification of our software-controlled HPLC and the peak parking procedure which will be used for any further separations.

John Percy for his programming skills used for the creation of the MoWeD deconvolution software which was initially used in the protein profiling experiments.

Roger Moore for creation of the winnow program which was used as an effective screen for MS/MS spectra.

This work was conducted under the Director's Research Discretionary Fund, Jet Propulsion Laboratory, California Institute of Technology