



# **ULTRA-SENSITIVE SENSOR BASED ON WHISPERING GALLERY MODE MICRO-RESONATORS**

Gregory Bearman, Vladimir Iltchenko and Lute Maleki

Jet Propulsion Laboratory  
California Institute of Technology  
Pasadena, CA

## Microsphere -- a low-loss photon trap, novel optical (micro)cavity

Whispering-gallery modes - closed circular waves under total internal reflection

(Term by J.W.S.Rayleigh, analogy to acoustic modes in the gallery of St Paul cathedral)

Sustained in any axisymmetric dielectric body with  $R \geq \lambda$

low material loss (transparent material, e.g fiber grade silica)

low bending loss ( $R \gg \lambda$ )

low scattering loss (TIR always under grazing incidence

+ molecular-size surface roughness)



Quality-factor $Q = \lambda/\Delta\lambda_{\text{RES}}$	-	up to $\sim 10^{10}$
---	---	----------------------

Photon lifetime $\tau = \lambda Q / 2\pi c$	-	up to $\sim 3\mu\text{s}$
---	---	---------------------------

(cavity ringdown time)


visible and near-infrared band: *Opt.Lett.* 21, p.453 (1996)

*Opt.Lett.* 23, p.247 (1998)



Visualization of WG mode field by residual scattering in silica microsphere, V.S.Ilchenko et al, *Opt.Comm.* 113, p.133(1994)

# Microsphere cavity modes sense external medium through evanescent wave





# Microspheres as sensors

Since evanescent wave samples immersion medium, losses will affect measured cavity  $Q$

Sensitivity can be enhanced by using absorption features to increase losses and further reduce  $Q$

As with immunoassays, use a fluorophore to indicate presence of analyte





A microsphere supported at the end of the glass rod on the right and optically coupled to the center fiber

# Measurement Technique


$Q < 10^8$

Since cavity bandwidth is larger than laser linewidth, we can scan the laser frequency and measure the  $Q$  directly


$Q > 10^8$

Use cavity ringdown to measure  $Q$


## Direct measurement of $\mu$ sphere Q









600  $\mu\text{m}$  sphere in 150 nM solution of TO-PRO-3 in  
Methanol  $\Delta=34$  MHz






600  $\mu\text{m}$  sphere in 500 nM solution of TO-PRO-3 in Methanol.  $\Delta=58$  MHz. Note the presence of other cavity modes that show up as the detuning range is  $\sim 1\text{GHz}$



Cavity ringdown measurement of microsphere cavity  $Q$ .  
 $Q = 1.2 \times 10^8$  vs  $1.1 \times 10^8$  from linewidth measurement



Measured concentrations of analyte with microsphere sensor