

Nanogeology of Metamorphic Magnetite using Three-Dimensional Atom Probe and Focused Ion Beam Secondary Ion Mass Spectroscopy

Kimberly R. Kuhlman, Pamela G. Conrad
California Institute of Technology, Jet Propulsion Laboratory

Richard L. Martens & Thomas F. Kelly,
Imago Scientific Instruments Corporation,
University of Wisconsin - Madison

Derren Dunn,
International Business Machines, Corp.

Neal D. Evans & Michael K. Miller
Oak Ridge National Laboratory

MEPAG/CAPTEM Workshop


August 15-16, 2001



<http://minileap.jpl.nasa.gov/>



California Institute of Technology



Objectives of this work

- A 3D near-atomic scale elemental map of a geological sample.
- Demonstrate the planetary science potential of the Local Electrode Atom Probe (LEAP).
 - Potential of LEAP analysis for non-conductive samples:
 - ➔ Terrestrial geology and geomicrobiology,
 - ➔ Apollo samples from the Moon,
 - ➔ Samples returned by the Stardust, Genesis and Mars missions.
- Demonstrate the potential of the Mini-LEAP for the *in-situ* analysis of planetary materials.
 - NASA is currently developing a prototype Mini-LEAP at JPL.





Schematic Illustration of APFIM Analysis




<http://minileap.jpl.nasa.gov/>

From Miller, et al. (1996) Atom Probe Field Ion Microscopy Oxford University Press



California Institute of Technology

Analytical Capability of 3DAP



Atom probe analysis in the vicinity of an interface in a pressure vessel steel. Enrichments of P, C, Mo and Mn are evident.



From Miller (2000) Atom Probe Tomography : Analysis at the Atomic Level, Plenum Publishing Corporation.

<http://minileap.jpl.nasa.gov/>



California Institute of Technology




Why Magnetite?

- A common mineral on Earth and Mars.
- One of the more conductive minerals.
 - Resistivity = 52×10^{-4} ohm-cm.
- This particular magnetite contains disk-shaped exsolutions approx. 40 nm in diameter, 1-3 nm thick and about 10^4 platelets/ μm^3 .
- EDS shows Mn and Al concentrated in these precipitates.
- Quantitative analysis has been limited by the thickness of this second phase.

Atom Probe is a technique that can potentially quantify the composition of these precipitates.



Recovery of a Magnetite Field Emission Tip



TEM of sample 031300E after resharpener with a Gatan ion mill with 5 keV Ar and gun tilt = 40°. a) Bright field image showing precipitates and bend contour. b) High magnification image of the tip and precipitates seen in FIM images. The line indicates the surface obtained after field evaporation.




<http://minileap.jpl.nasa.gov/>




California Institute of Technology

Field Ion Micrographs of Magnetite LP204-1



16 kV



16.5 kV

Note the clarity of the $\{111\}$ poles and the precipitates running directly through the pole at right of the images.




<http://minileap.jpl.nasa.gov/>




California Institute of Technology

Six-fold Symmetry of the $\{111\}$ Planes



16 kV



16.5 kV




<http://minileap.jpl.nasa.gov/>




California Institute of Technology

Post Evaporation TEM Analysis of Magnetite LP204-1



Disk-shaped precipitates
normal to the beam.



Precipitates seen in FIM
images. The location of
this surface is indicated on
previous slide.

Electron diffraction pattern
indicates that this view is
normal to the [100].




<http://minileap.jpl.nasa.gov/>



California Institute of Technology



First Mass Spectra from a Geological Material using Atom Probe




<http://minileap.jpl.nasa.gov/>




California Institute of Technology



Preliminary FIB-SIMS Imaging of Magnetite



Secondary Ion Image



Aluminum Map



<http://minileap.jpl.nasa.gov/>



California Institute of Technology




FIB-SIMS Image of a Cross-Sectioned Spore



<http://minileap.jpl.nasa.gov/>



California Institute of Technology



Conclusion

- First high quality FIM images of a bulk mineral.
- Preliminary APFIM analysis of magnetite:
 - Mass spectra have been obtained showing isotopes of Fe and O as well as FeO.
 - Singly and doubly charged ions are observed.
 - A preliminary mass spectrum taken from a precipitate contains both Mn and Al, as expected from previous analysis of magnetite LP204-1 by Sitzman, et al.
- Traditionally non-conductive samples **CAN** be imaged and analyzed using APFIM.
- FIB-SIMS can provide high-resolution chemical imaging of minerals.





Acknowledgements

➤ *Research sponsored by:*

- Cross Enterprise Technology Development Program, U.S. National Aeronautics and Space Administration.
- Director's Research and Development Fund, Jet Propulsion Laboratory.
- Oak Ridge National Laboratory SHaRE User Facility and SHaRE Program.
 - ➔ *Division of Materials Sciences and Engineering, U.S. Department of Energy.*
 - ➔ *SHaRE Program with Oak Ridge Associated Universities.*

➤ *Special thanks to:*

- Professor John Valley of the University of Wisconsin for providing the samples of magnetite used in this work.
- William Carmichael and Glenn Boda of the Electron Microscopy Program at the Madison Area Technical College in Madison, WI for their assistance and for providing access to the focused ion beam.

